



Equipment

- Aluminium construction
- Handle height adjustable 86 – 97 mm
- Incl. stick holder and big fabric bag

Data

Max. capacity:	136 kg
Total weight:	5,8 kg
Seat height:	54 cm
Seating area:	45 x 25 cm
Total length:	76,5 cm
Total width:	64 cm
Total height:	86 – 97 cm
Color:	Platin, Midnight Blue

Price

259,90 Euro (RRP)

Light weight rollator
CAPERO

With the CAPERO you will explore your everyday life newly and free. The Aluminum rollator weighs only 5,8 kg, including some extras which are serial equipment.

Capero is equipped with a stick holder and a bag for personal belongings or for small purchases.

If you want to have a little break, the nylon fabric and a stabile backrest belt allow comfortable seating. The height adjustable push handles enable Capero to adjust to your individual body measures.

Transport is also very easy with this light weight. With one pull, it folds together and additionally, with one click, wheels can be removed, if needed. Like this CAPERO will fit into any car without problem.

Characteristics:

- weather-proof bag (max. capacity 5 kg)
- Total weight only 5,8 kg
- removable wheels