



## Light weight rollator **RL-SMART**

The light weight rollator out of aluminium weighs 6.7 kg only. A specially broad and soft seat enables you to have comfortable rest periods.

Reflectors on front wheels provide for safety in the dark. Due to push handles adjustable in height the rollator can be used by tall and small people.

Delivery contains a stick holder as well as a large removable fabric bag (max. capacity 5 kg) as an additional extra. The new light weight rollator visually shines – in platinum or blue metallic colour- and with a puristic design, which is underlined by inner tube brake cables.

### Equipment

Aluminium construction  
Handle height adjustable 82 – 94,5 cm

Stick holder

Big fabric bag

### Data

Max. capacity:	136 kg
Total weight:	6,7 kg
Seat height:	54 cm
Seating area:	45 x 25 cm
Total length:	65,5 cm
Total width:	63 cm
Total height:	82 – 94,5 mm
Color:	Platin   Blue metallic

### Price

**195,00 Euro** (RRP)

### Characteristics:

- Specially broad and soft seat
- Light weight due to aluminium construction
- Puristic design