

**TOPRO
OLYMPOS**

TOPRO



THE OFFROADER!
A tough and practical rollator for an active life

www.topro.co.uk

**TOPRO
OLYMPUS**





Designed especially for outdoor use, but can also be used inside

- The tyres are made for comfortable driving on rough foundation. For example cobblestones or dirt roads.
- Spacious net basket with a slipcover.
- Made from heavy duty materials which can tolerate all kinds of weather.
- Easy to clean.
- Tested and approved according to EN ISO 11199



Easy to unfold and a height adjustment with memory function

- Unique memory function that sets the handles back in the right height.



Safe to use with excellent brakes, kerb lifter and edge guard.

- Safe to use and practical when crossing thresholds, kerbs etc.
- Easy to operate brakes.
- Anatomic handles.
- A variety of accessories that give extra safety in different situations.



Folds easy together and requires little space

- The rollator folds together by pulling the strap on the seat.
- Stands steady on its wheels when folded.
- Takes little space when being stored.
- Is practical and easy to transport.
- Excellent for users who travel by car, bus, train or airplane.



Size	M (medium)		S (small)	
	KG	LBS	KG	LBS
Colour	Dark grey		Silver	
Maximum weight of user	150	330	150	330
Product weight	8,6	19	8,0	17 lbs 6 oz
Weight capacity basket	15	33	15	33
	CM	INCHES	CM	INCHES
Recommended user height	150 - 195	4ft 11" - 6ft 4 3/8"	135 - 170	4ft 5 1/8" - 5ft 7"
Height of handles	76 - 95	30" - 37 3/8"	70 - 89	27 1/2" - 35"
Maximum length	76	30"	72	28 3/8"
Maximum width	70,5	27 3/4"	66	26
Width between handles	50	19 3/4"	46	18 1/8"
Seat height	60	23 5/8"	53	20 7/8"
Turning diameter	94	37"	86,5	34"
Diameter front wheels	27,5	10 7/8"	22,5	8 7/8"
Diameter rear wheels	22,5	8 7/8"	20	7 7/8"
Width on front/rear wheels	4,2	1 5/8"	4,2	1 5/8"
Grip distance driving brake	6,5	2 5/8"	6,5	2 5/8"
Height when folded up	86	33 7/8"	78	30 3/4"
Width when folded up	27	10 5/8"	27	10 5/8"
Length when folded up	76	30	72	28 3/8"
TOPRO ART. No.	814300		814305	

Accessories for a variety of uses

To be ordered separately.

Product **TOPRO ART. No.**



Adjustable back support with padding, for rollators from 2009 **814726**



Crutch holder **814040**



Bottle holder **814043**



Rear bag with zipper **814046**



Transport bag **814042**



New!
LED light with integrated fall alarm (pair) **814639**



Lamp **814608**



Nameplate **814024**



Lock **814023**



Ice spurs for shoes
Size M (shoe size 36-40) **814029**
Size L (shoe size 41-46) **814035**



Umbrella with attachment **814068**



Basket for oxygen bottle (only compressed, not liquid) **814009**



Guiding handle for the visually impaired **814014**



One-handed brake **814026**



New!
Drag brake (pair) **814032**