

TOPRO
TROJA Classic

TOPRO



THE ORIGINAL!

A supreme test winner on the market since 2002

The acknowledged German test institute, Stiftung Warentest, tested 14 rollators in edition 9/2005. This is what they said about the Norwegian manufactured test winner, TOPRO TROJA:

- Best operating comfort.
- Easy to fold.
- Good brake function.
- Light weight and easy to manoeuvre.
- Requires little space when folded.
- The test team's favorite.

www.topro.co.uk



**TOPRO TROJA ROLLATORS
are easy to fold, requiring little space when stored**

- Pull the cord on the seat and the rollator folds together horizontally.
- Push the handlebars down into the sidebars.
- Unique memory function that sets the handles back to the right height.
- Stands securely on its wheels when folded.
- Takes little space and is very practical to transport.
- Excellent for users who travel by car, train, bus or airplane.

An elegant and practical light weight rollator

- Plastic coated aluminium tube and plastic couplings.
- A total of 12 ball bearings in the wheels and forks provide excellent driving comfort.
- Can be hosed and cleaned.
- All wear and tear parts can be replaced and brakes can be adjusted.
- Tested and approved according to EN ISO 11199

Excellent brakes, kerb lifter and edge guard

- Easy to negotiate thresholds and kerbs etc.
- Safe to use.

Three sizes!
Colour

Maximum weight of user
Product weight

Recommended user height
Height of handles Maximum
length Maximum width

Width between handles

Seat height
Turning diameter

Diameter front/rear wheels
Width front/rear wheels Grip

distance driving brake
Height when folded up

Width when folded up
Length when folded up

Handles
Area of use

TOPRO ART. No.


M (medium)

Dark grey	
KG	LBS
165	364
7,4	16,3
CM	INCHES
150-200	4ft 11"-6ft 4 3/8"
78-100	30 3/4"-39 3/8"
65	25 5/8"
61	24
47	18 1/2"
62	24 3/8"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
80	31 1/2"
24	9 1/2"
66	26

S (small)

Silver	
KG	LBS
125	274
7,1	15,6
CM	INCHES
135-170	4ft 5 1/8"-5ft 7"
67-86	26 3/8"-33 1/8"
65	25 5/8"
61	24
47	18 1/2"
54	21 1/4"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
74	29 1/8"
24	9 1/2"
66	26

X (extra small)

Silver	
KG	LBS
80	176
6,6	14,2
CM	INCHES
125-160	4ft 1 1/2"-5ft 3"
63,5-79	25"-31"
65	25 5/8"
51,5	20 1/4"
37	14 1/2"
52	20 1/2"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
65	25 5/8"
25	9 3/8"
66	26

Anatomic

Anatomic

Anatomic

Designed and approved for both indoor and outdoor use

814750
814770
814740

Also available complete with forearm supports. See TOPRO Walkers brochure.

Accessories for a variety of uses

To be ordered separately.

Product **TOPRO ART. No.**

Shopping bag with cover and shoulder strap **X** **814140**



Ice spurs for shoes
Size M (shoe size 36-40) **814029**
Size L (shoe size 41-46) **814035**



Adjustable back support with padding, for rollators from 2009 **814726**



Rear bag with zipper **X** **814046**



Umbrella (not for sale in Germany) **814735**
157207



Back Support for TROJA Walker **814765**



Transport bag **814042**



Lock **814023**



Tray **X** **814728**



Flower stickers, set with 4 stickers **814141**



Guiding handle for the visually impaired **814014**



Anti slip mat for tray **X** **814059**



Reflective Panels **814624**



Basket for oxygen bottle (only compressed, not liquid) **X** **814009**



Crutch holder **814025**



Lamp **814608**



One-handed brake **814026**
157197



Bottle holder **814043**



Nameplate **814024**



Forearm support (pair) for TROJA M **814829**
for TROJA S and TROJA X **814845**