



Elegant and practical light weight walkers

- Plastic coated aluminium tube and plastic couplings.
- A total of 12 ball bearings in the wheels and forks provide excellent driving comfort.
- Can be hosed and cleaned.
- All wear and tear parts can be replaced and brakes can be adjusted.
- Tested and approved according to EN ISO 11199-3

Easy to fold, requiring little space when stored

- Pull the cord on the seat and the walker folds together horizontally.
- Push the handle bars down into the side bars.
- Unique memory function that sets the handles back to the right height.
- Stands securely on its wheels when folded.
- Takes little space and is very practical to transport.
- Excellent for users who travel by car, train, bus or airplane.

Excellent brakes, kerb lifter and edge guard

- Easy to negotiate thresholds and kerbs etc.
- Safe to use.



Three sizes !

M (medium)

S (small)

X (extra small)

Colour	M (medium)		S (small)		X (extra small)	
	KG	LBS	KG	LBS	KG	LBS
Maximum weight of user	135	298	125	274	80 kg	176
Product weight	10,1	22	9,8	21,8	9,3 kg	20,6
	CM		INCHES		CM	
Recommended user height	150-185	4ft 11"-6ft 1"	130-170	5ft 1 1/4"-5ft 7"	130-160	5ft 1 1/4"-5ft
Height of handles	98-115	38 5/8"-45 1/4"	84-100	33"-9 3/8"	84-96	33"-37 3/4"
Maximum length	65	25 5/8"	65	25 5/8"	65	25 5/8"
Maximum width	61	24	61	24	52,5	20 1/2"
Width between handles	45	17 3/4"	45	17 3/4"	35	13 3/4"
Seat height	62	24 3/8"	54	21 1/4"	52	20 1/2"
Turning diameter	84	33 1/8"	84	33 1/8"	84	33 1/8"
Diameter front/rear wheels	20	7 7/8"	20	7 7/8"	20	7 7/8"
Width front/rear wheels	3,6	1 3/8"	3,6	1 3/8"	3,6	1 3/8"
Diameter handles (round)	3,5	1 3/8"	3,5	1 3/8"	3,5	1 3/8"
Grip distance driving brake	6,5	2 1/2"	6,5	2 1/2"	6,5	2 1/2"
Height when folded up	108	42 1/2"	88	34 5/8"	88	34 5/8"
Width when folded up	31	12 1/4"	31	12 1/4"	31	12 1/4"
Length when folded up	66	26	66	26	66	26

Designed and approved to be used indoors on plain foundation

TOPRO ART. No.

814749

814742

814741

Imperial measurements are approximate.

Accessories for a variety of uses

To be ordered separately.

Product **TOPRO ART. No.**



Back Support **814765**



Bottle holder **814043**



Nameplate **814024**



Tray X **814728**



Shopping bag with cover and shoulder strap X **814140**



Lock **814023**



Anti slip mat for tray X **814059**




Rear bag with zipper X **814046**



Guiding handle for the visually impaired **814014**



Crutch holder **814025**



Flower stickers, set with 4 stickers **814141**



Basket for oxygen bottle (only compressed, not liquid) X **814009**